

Lovibond® Colour Measurement

Tintometer® Group



TRA 520 / TRA 500

- Portable docking station for use with the Lovibond® TR 520 / 500 series of spectrophotometers
- Provides a light controlled environment and sample optical cells for the colour measurement of liquids, gels, powders and other 'awkward' objects
- No ambient light - tight-fitting integral lid
- Supports a wide variety of samples with cell sizes of 10, 20 and 30mm
- Supplied with a white background plate to ensure consistency of readings
- Light weight and portable, yet very sturdy and robust
- Quick and easy to calibrate the instrument within the adapter
- Available as a package - TR 500 (TRA 500) and TR 520 (TRA 520) or as an upgrade from the TR 500 / 520 (serial number required)

Market sectors

Coatings

Food and Beverage

Household Products


Medical and Clinical

Plastics




TRA 520 / TRA 500 Innovative Product Design

The TRA 500 / 520 docking station is precision engineered to house the TR500 / TR 520 spectrophotometer at the optimum ergonomic position for touch screen operation and sample colour measurement.



A cover slots securely into place above the glass cell to ensure no ambient light can affect the reading



10, 20 and 30mm optical glass cells are supplied with the unit to securely position and hold a variety of liquid, gel and powder samples

A white background plate is included to ensure consistency of colour measurement

Order No	Product	Order No	Product
403225	Lovibond® TRA 520	605910	W586.OG.30mm Cell
403220	Lovibond® TRA 500 (8mm aperture)	605895	W586.OG.10mm Cell (Pack of 5)
605890	W586.OG.10mm Cell	605905	W586.OG.20mm Cell (Pack of 5)
605900	W586.OG.20mm Cell	605915	W586.OG.30mm Cell (Pack of 5)

 www.lovibond.com

 sales@tintometer.com

 [Download form](#)

Sphera Analytics GmbH

Keltenring 1-3, 82041 Oberhaching, Germany

Email: info@sphera-analytics.com

Web: <https://sphera-analytics.com>

Lovibond® and Tintometer® are Registered Trade Marks of The Tintometer® Group. All translations and transliterations of LOVIBOND® & TINTOMETER® are asserted as Trade Marks of The Tintometer Limited. Registered Office: Lovibond House, UK. Registered in England No. 45024.

Errors and Omissions Excepted – Content Subject to alterations without notice.

The Bluetooth® mark, figure mark and combination mark are trademarks owned by the Bluetooth SIG and any use by the Tintometer® Group is under license. Please note, Bluetooth® is not available for all regions. Please contact the local Lovibond® representative for details.